


# CT Dose Profiler



## A Pioneer at Measuring CT Dose

The RTI CT Dose Profiler has taken the CT quality assurance to the next level. Because of its revolutionary design it has transformed the CTDI measurement from being inaccurate due to underestimation of the dose for wide beams to be more exact. It also has the ability to further analyze the result.

All in one shot.

## Specifications

### CT Dose Profiler

**Supporting Barracuda electrometer:**

All electrometer modules

**Max sensitivity variation (0°-360°):**

less than  $\pm 5\%$

**Typical calibration factor:** 0.3 mGy/nC

**Material:** Al and PMMA

**Connector:** Triaxial LEMO

**Diameter:** 12.5 mm

**Detector width:** 0.3 mm

**Length:** 165 mm + 45 mm

**Trig modes:** Timed, After Exposure and Continuous

**Max scanning time:** 160 sec

*Innovative X-ray QA Solutions...of Course!*




# A Small Piece of Revolution

## No Limitations

Today computed tomography (CT) contributes up to 70% of the total dose given to patients during X-ray examinations. The rapid advancements in CT technology are placing new demands on the methods and equipment that are used for quality assurance. The wide beam widths found in CT scanners with multiple beam apertures make it impossible to use existing CT ionization chambers to measure the total dose given to the patient. Using a standard 10 cm CT ionization chamber may result in inaccurate measurements due to underestimation of the dose profile for wide beams. The CT Dose Profiler Probe is a new type of CT detector that does not have this limitation. It can be used with the Barracuda or the Piranha and a PC running the CT Dose Profile Analyzer Software.

The CT Dose Profiler is based on solid-state technology. It is robust and it fits into existing standard phantoms used for CTDI measurements. The sensitive part of the detector chip is very thin (width 0.3 mm). Thanks to the small width, the detector is completely irradiated when the table is moving and the CT scans over the probe. The dose is measured in every point of the X-ray beam and the total dose profile is acquired regardless of beam width. This makes it possible to measure without the drawbacks of traditional CT probes.



## Measure and Analyze

Following parameters are achieved from a single exposure:

- $CTDI_{100}$
- $CTDI_w$
- $CTDI_{vol}$
- DLP
- Performance of the AEC
- Geometric efficiency (%)
- FWHM (Full width at half maximum of the dose profile.)
- Scatter Index (An index for checking how much of the total dose that is not measured by a standard CT ionization chamber.)
- CT dose profile
- Point Dose



### CT Phantom

Phantoms for dose measurements on CT scanners. Body and head Phantom with insert comes in a hard case with built-in trolley.



# Scandinavian Quality

RTI Electronics was founded in 1981 when several curious and enterprising students met at Chalmers University of Technology in Gothenburg, Sweden. They saw their vision grow into the beginning of RTI products – today world leading in X-ray QA and Service instrumentation.

There are many reasons why RTI Electronics has become a market leader. Besides fulfilling the highest user demands, products from RTI Electronics are known for cutting edge innovation. Other reasons include our engagement, our expertise accumulated over more than a quarter of a century, and our commitment to doing it right.

We are convinced that You will be satisfied with Your choice of product, and we would like to continue to grow – together with You.

© Copyright 2009 RTI Electronics AB – CT Dose Profiler BAR PIR 201005



#### World Headquarters

RTI Electronics AB  
Flöjelbergsgatan 8 C  
SE-431 37 Mölndal  
SWEDEN

Phone: + 46 31 746 36 00  
Fax: + 46 31 27 05 73  
E-mail: [sales@rti.se](mailto:sales@rti.se)  
[www.rti.se](http://www.rti.se)

#### US Office

RTI Electronics, Inc.  
1275 Bloomfield Avenue  
Building 5, Unit 29A  
Fairfield, NJ 07004  
USA

Phone: 800-222-7537  
Phone: 1-973-439-0242  
Fax: 1-973-439-0248  
E-mail: [sales@rtielectronics.com](mailto:sales@rtielectronics.com)  
[www.rtielectronics.com](http://www.rtielectronics.com)