

Field Position Analyzer

The Piranha Panoramic Holder with straps is specially designed for attaching the Piranha to dental panoramic equipment. With this holder you'll get an easy and quick positioning of the Piranha without any hassle with tape.



World Headquarters

RTI Electronics AB
Flöjelbergsgatan 8 C
SE-431 37 Mölndal
SWEDEN

Phone: + 46 31 746 36 00
Fax: + 46 31 27 05 73
E-mail: sales@rti.se
www.rti.se

US Office

RTI Electronics, Inc.
33 Jacksonville Road, Bldg. 1
Towaco, NJ 07082
USA

Phone: 1-800-222-7537
Fax: 1-973-439-0248
E-mail: sales@rtielectronics.com
www.rtielectronics.com

Field Position Analyzer that Saves Time and Money

The Piranha Panoramic Holder with straps is specially designed for attaching the Piranha to dental panoramic equipment. With this holder you'll get an easy and quick positioning of the Piranha without any hassle with tape. The fluorescent markings with scales guides you to correct positioning. The holder can also with advantage be used for other purposes where the Piranha has to be attached to a surface, e.g. Cone Beam CT units.

Note that on panoramic systems with sliding cassette, the Piranha holder with straps is not possible to use. In these cases we recommend the *Piranha Panoramic Holder Magnetic* instead.

Piranha Measurement on Panoramic Systems

The Piranha is ideal for measurement on panoramic (OPG) systems and all other vertical positioning of the Piranha. The small detector area and the "Piranha position check", a feature that makes it possible to verify the detector position, ensures more accurate and more reproducible results than with other devices.

Since the Piranha works wireless, you do not have to worry about any cable during the rotation of the OPG. The Piranha Panoramic holder with straps makes it easy to attach the Piranha on the panoramic system and quickly perform the measurement. Of course, the Piranha is perfectly suited for measurements on standard intra oral systems as well.



How to Use the Panoramic Holder with Straps for Piranha video.

To Perform a Measurement

- Place the Piranha inside the "Snap-in" Piranha Holder.
- Attach the suction cup to help secure the Panoramic Holder and strap.
- Wrap the Multi-adjustable strap around the machine housing and lock the clips. Tighten if necessary.
- Make sure the Piranha is in position by sliding it to the right or to the left. The fluorescent screen helps guide you to the desired position.

Field Position Analyzer



The Piranha Panoramic Holder Straps consists of the following parts:

1. Easy "Snap-In" Piranha Holder
2. Fixation Suction Cup
3. Multi-adjustable Strap with locking clips