


<p>Modality</p> <p>Measurements</p>	Intraoral	Panoramic	Cone Beam CT
Quick Check	<p>Cobia Dental</p> 	<p>Piranha 205</p> 	<p>Cobia Dental</p> 
Quick Check & Service	<p>Piranha 205</p> 	<p>Piranha 205</p> 	<p>Piranha 205</p> 
Complete X-ray Quality Assurance	<p>Piranha Premium Kit</p> 	<p>Piranha Premium Kit</p> 	<p>Piranha Premium Kit</p> 
Suggested Accessories	 <p>Panoramic Holder provides an easy positioning for measurement on Panoramic systems.</p>		 <p>The RTI Mover makes it possible to measure dose profiles with the RTI CT Dose Profiler and the Piranha.</p>